






Feature	EpiPen®	Anapen®
Adrenaline dose	Single pre-measured	Single pre-measured
150 microgram device 7.5-20kg (aged around 1-5 years) Green label	 <p>EpiPen® Junior (150 microgram)</p>	 <p>Anapen® Jr (150 microgram)</p>
300 microgram device Children 20kg and over (aged around 5 years and over) and adults Yellow label	 <p>EpiPen® (300 microgram)</p>	 <p>Anapen® 300 (300 microgram)</p>
500 microgram device Children 50kg and over (aged around 12 years and over) and adults Magenta (purple) label	Not available	 <p>Anapen® 500 (500 microgram)</p>
Viewing window to check adrenaline for discolouration or precipitate	Yes	Yes
Availability	Available over-the-counter at full price 2 devices on Pharmaceutical Benefits Scheme (PBS) authority prescription	Available over-the-counter at full price 2 devices on Pharmaceutical Benefits Scheme (PBS) authority prescription
Activation of device	Press firmly against outer mid-thigh	Press red button when device on outer mid-thigh
Safety	Blue safety release Orange needle end automatically extends over needle after use	Grey safety cap Black needle shield can be replaced over needle after use
Trainer devices	Available from Allergy & Anaphylaxis Australia - https://allergyfacts.org.au/shop/adrenaline-injector-trainers-and-resources	Available from Allergy & Anaphylaxis Australia - https://allergyfacts.org.au/shop/adrenaline-injector-trainers-and-resources
Expiry reminder service	MyEpiPen https://www.myepipen.com.au/	Optional patient sms/email Anapen® expiry reminder available at participating pharmacies https://anapen.pharmaprograms.com.au/locator/